

## Product Datasheet

### Human CCL1 ELISA Kit (orb405664)

<b>Description</b>	Human C-C motif chemokine 1(CCL1) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6 pg/ml - 1000 pg/ml
<b>Target</b>	CCL1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P22362</a>
<b>Sensitivity</b>	3.9 pg/ml
<b>Expiration Date</b>	Please enquire.