

Product Datasheet

Human CCR2 ELISA Kit (orb405667)

Description	Human CC-Chemokine Receptor 2,CCR2 ELISA Kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	CCR2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid (CSF)
Uniprot ID	P41597
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.