

Product Datasheet

Human CXCR4 ELISA Kit (orb405703)

| | |
|--------------------------|--|
| Description | Human CXC-chemokine receptor 4(CXCR4)ELISA Kit |
| Reactivity | Human |
| Range | 1.56 pg/mL-100 pg/mL |
| Target | CXCR4 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P61073 |
| Sensitivity | 0.39 pg/mL |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)