

## Product Datasheet

### Human CXCR4 ELISA Kit (orb405703)

<b>Description</b>	Human CXC-chemokine receptor 4(CXCR4)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6 pg/mL-1000 pg/mL
<b>Target</b>	CXCR4
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P61073</b>
<b>Sensitivity</b>	3.9 pg/mL
<b>Expiration Date</b>	Please enquire.