

Product Datasheet

Human FLT3LG ELISA Kit (orb405734)

Description	Human FMS-like tyrosine kinase 3 ligand,Flt-3L ELISA Kit
Reactivity	Human
Range	7.8 pg/mL-500 pg/mL
Target	FLT3LG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	P49771
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.