

## Product Datasheet

### Human ICOSLG ELISA Kit (orb405754)

<b>Description</b>	Human ICOS ligand(ICOSLG) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	93.75 pg/mL-6000 pg/mL
<b>Target</b>	ICOSLG
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">O75144</a>
<b>Sensitivity</b>	23.4 pg/mL
<b>Expiration Date</b>	Please enquire.