

Product Datasheet

Human LDHB ELISA Kit (orb405770)

Description	Human L-lactate dehydrogenase B chain(LDHB) ELISA kit
Reactivity	Human
Range	23.5 mU/mL-1500 mU/mL
Target	LDHB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid (CSF), urine.
Uniprot ID	P07195
Sensitivity	5.8 mU/mL
Expiration Date	Please enquire.