

## Product Datasheet

### Human S100A12 ELISA Kit (orb405820)

**Description** Human S100 calcium binding protein A12/Calgranulin-C(S100A12) ELISA Kit

**Reactivity** Human

**Range** 0.156 ng/ml-10 ng/ml

**Target** S100A12

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [P80511](#)

**Sensitivity** 0.039ng/ml.

**Expiration Date** Please enquire.