

Product Datasheet

Human S100A9 ELISA Kit (orb405824)

| | |
|--------------------------|--|
| Description | human S100 calcium binding protein A9/calgranulin B,S100A9 ELISA Kit |
| Reactivity | Human |
| Range | 4.69 ng/mL-300 ng/mL |
| Target | S100A9 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P06702 |
| Sensitivity | 1.17 ng/mL |
| Expiration Date | Please enquire. |