

## Product Datasheet

### Human SMOC2 ELISA Kit (orb405841)

<b>Description</b>	Human SPARC-related modular calcium-binding protein 2(SMOC2) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Target</b>	SMOC2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9H3U7</a>
<b>Sensitivity</b>	0.039 ng/mL
<b>Expiration Date</b>	Please enquire.