

Product Datasheet

Human TIMD4 ELISA Kit (orb405869)

Description	Human T-cell immunoglobulin and mucin domain-containing protein 4(TIMD4) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	TIMD4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q96H15
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.