

## Product Datasheet

### Human alpha 1 Antitrypsin ELISA Kit (orb405880)

<b>Description</b>	Human $\alpha$ 1-Antitrypsin, $\alpha$ 1-AT ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78 nIU/mL-5000 nIU/mL
<b>Target</b>	SERPINA1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Uniprot ID</b>	<b>P01009</b>
<b>Sensitivity</b>	19.5 nIU/mL
<b>Expiration Date</b>	Please enquire.