

Product Datasheet

Human SERPINA3 ELISA Kit (orb405881)

Description	Human Alpha1 Antichymotrypsin,AACT ELISA Kit
Reactivity	Human
Range	1.25 ng/ml-80 ng/ml.
Target	SERPINA3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P01011
Sensitivity	0.31 ng/ml
Expiration Date	Please enquire.