

## Product Datasheet

### Human $\alpha$ -GST ELISA Kit (orb405895)

<b>Description</b>	Human $\alpha$ -glutathione S-transferases, $\alpha$ -GST ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Target</b>	$\alpha$ -GST
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P08263</b>
<b>Sensitivity</b>	0.039 ng/mL
<b>Expiration Date</b>	Please enquire.