

## Product Datasheet

### Human bCTx ELISA Kit (orb405927)

<b>Description</b>	Human beta-crosslaps (bCTx) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	125 pg/mL-2000 pg/mL
<b>Target</b>	bCTx
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Sensitivity</b>	125 pg/mL
<b>Expiration Date</b>	Please enquire.