

## Product Datasheet

### Human FGF10 ELISA Kit (orb406114)

<b>Description</b>	Human Fibroblast Growth Factor 10 (FGF10) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.12 pg/mL-200 pg/mL
<b>Target</b>	FGF10
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">O15520</a>
<b>Sensitivity</b>	0.78 pg/mL
<b>Expiration Date</b>	Please enquire.