

Product Datasheet

Human CD274 ELISA Kit (orb406126)

Description	Human Programmed Death Ligand-1(PD-L1/CD274) ELISA Kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	CD274
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9NZQ7
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.