

## Product Datasheet

# Human anti-Survivin(Surv) antibody ELISA Kit (orb406149)

<b>Description</b>	Human anti-Survivin(Surv) antibody ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625 ng/mL-40 ng/mL
<b>Target</b>	Surv Ab
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Sensitivity</b>	0.156 ng/mL
<b>Expiration Date</b>	Please enquire.