

## Product Datasheet

### Human CCL13 ELISA Kit (orb406157)

<b>Description</b>	Human monocyte chemotactic protein 4, MCP-4 ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	39 pg/ml-2500 pg/ml.
<b>Target</b>	CCL13
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates.
<b>Uniprot ID</b>	<a href="#">Q99616</a>
<b>Sensitivity</b>	9.75 pg/ml.
<b>Expiration Date</b>	Please enquire.