

Product Datasheet

Human PIN1 ELISA Kit (orb406168)

Description	Human protein(peptidylprolyl cis/trans isomerase)NIMA-interacting 1,PIN1 ELISA Kit
Reactivity	Human
Range	1.56 pg/mL-100 pg/mL
Target	PIN1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, tissue homogenates
Uniprot ID	Q13526
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.