

Product Datasheet

Human Protein S ELISA Kit (orb406171)

| | |
|--------------------------|-----------------------------------|
| Description | Human Protein S ELISA Kit |
| Reactivity | Human |
| Range | 3.12 ng/mL-200 ng/mL |
| Target | PROS1 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | P07225 |
| Sensitivity | 0.78 ng/mL |
| Expiration Date | Please enquire. |