

## Product Datasheet

### Human LRP4 ELISA Kit (orb406191)

<b>Description</b>	Human low-density lipoprotein-receptor-related protein 4(LRP-4)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5 pg/mL-4000 pg/mL
<b>Target</b>	LRP4
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>075096</b>
<b>Sensitivity</b>	15.6 pg/mL
<b>Expiration Date</b>	Please enquire.