

Product Datasheet

Human HIF-1 α ELISA Kit (orb406195)

Description	Human hypoxia-inducible factor 1 α (HIF-1 α) ELISA kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	HIF-1 α
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q16665
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.