

Product Datasheet

Human Versican ELISA Kit (orb406218)

| | |
|--------------------------|---|
| Description | human versican/PG-M/PG-350 ELISA kit |
| Reactivity | Human |
| Range | 0.04 ng/ml-30 ng/ml. |
| Target | VCAN |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates. |
| Uniprot ID | P13611 |
| Sensitivity | 0.04 ng/ml. |
| Expiration Date | Please enquire. |