

## Product Datasheet

# Human ReverseTri-iodothyronine ELISA Kit (orb406236)

<b>Description</b>	Human ReverseTri-iodothyronine,rT3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	75 pg/mL-1200 pg/mL
<b>Target</b>	rT3
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Sensitivity</b>	75 pg/mL
<b>Expiration Date</b>	Please enquire.