

## Product Datasheet

### Human SFTPD ELISA Kit (orb406246)

<b>Description</b>	Human Pulmonary surfactant-associated protein D,SP-D ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.25 ng/mL-80 ng/mL
<b>Target</b>	SFTPD
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, bronchoalveolar lavage fluid, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P35247</a>
<b>Sensitivity</b>	0.31 ng/mL
<b>Expiration Date</b>	Please enquire.