

## Product Datasheet

### Human Galanin ELISA Kit (orb406292)

<b>Description</b>	Human galanin,GAL ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5 pg/ml-4000 pg/ml.
<b>Target</b>	GAL
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates.
<b>Uniprot ID</b>	<b>P22466</b>
<b>Sensitivity</b>	15.6 pg/ml
<b>Expiration Date</b>	Please enquire.