

## Product Datasheet

### Human MASP1 ELISA Kit (orb406293)

<b>Description</b>	Human mannan-binding lectin serine peptidase 1 (C4/C2 activating component of Ra-reactive factor) (M...
<b>Reactivity</b>	Human
<b>Range</b>	3.12 ng/mL-200 ng/mL
<b>Target</b>	MASP1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P48740</a>
<b>Sensitivity</b>	0.78 ng/mL
<b>Expiration Date</b>	Please enquire.