

Product Datasheet

Human BMP2 ELISA Kit (orb406351)

Description	Human Bone morphogenetic protein 2,BMP-2 ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	BMP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, urine, cell culture supernates, tissue homogenates
Uniprot ID	P12643
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.