

## Product Datasheet

### Human PRDX6 ELISA Kit (orb406391)

<b>Description</b>	Human Peroxiredoxin-6(PRDX6) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	12.5 pg/mL-800 pg/mL
<b>Target</b>	PRDX6
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P30041</b>
<b>Sensitivity</b>	3.12 pg/mL
<b>Expiration Date</b>	Please enquire.