

Product Datasheet

Human KLK5 ELISA Kit (orb406480)

Description	Human Kallikrein-5(KLK5) ELISA kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	KLK5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, saliva
Uniprot ID	Q9Y337
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.