

## Product Datasheet

### Human HAVCR2 ELISA Kit (orb406491)

<b>Description</b>	Human Hepatitis A virus cellular receptor 2(HAVCR2) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	47 pg/mL-3000 pg/mL
<b>Target</b>	HAVCR2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q8TDQ0</a>
<b>Sensitivity</b>	11.75 pg/mL
<b>Expiration Date</b>	Please enquire.