

## Product Datasheet

### Human PTHLH ELISA Kit (orb406495)

**Description** Human Parathyroid Hormone Related Protein,PTHrP ELISA Kit

**Reactivity** Human

**Range** 31.25 pg/mL-2000 pg/mL

**Target** PTHLH

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [P12272](#)

**Sensitivity** 7.81 pg/mL

**Expiration Date** Please enquire.