

## Product Datasheet

### Human CTGF ELISA Kit (orb406522)

**Description** Human connective tissue growth factor (CTGF) ELISA kit

**Reactivity** Human

**Range** 3.75 pg/mL-120 pg/mL

**Target** CTGF

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [P29279](#)

**Sensitivity** 1.875 pg/mL

**Expiration Date** Please enquire.