

## Product Datasheet

### Human MIP-2 ELISA Kit (orb406547)

<b>Description</b>	Human macrophage inflammatory protein 2, MIP-2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.12 pg/mL-200 pg/mL
<b>Target</b>	MIP-2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, urine
<b>Uniprot ID</b>	<b>P19875</b>
<b>Sensitivity</b>	7.81 pg/mL
<b>Expiration Date</b>	Please enquire.