

Product Datasheet

Human AMHR2 ELISA Kit (orb406634)

Description Human Anti-Muellerian hormone type-2 receptor(AMHR2) ELISA kit

Reactivity Human

Range 0.32 ng/mL-12 ng/mL

Target AMHR2

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates, cell lysates

Uniprot ID [Q16671](#)

Sensitivity 8 pg/mL

Expiration Date Please enquire.