

## Product Datasheet

# Human myelin oligodendrocyte glycoprotein antibody (IgG) ELISA Kit (orb406658)

<b>Description</b>	Human myelin oligodendrocyte glycoprotein antibody (IgG) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	Request Information
<b>Target</b>	MOG Ab (IgG)
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum
<b>Sensitivity</b>	Request Information
<b>Expiration Date</b>	Please enquire.