

## Product Datasheet

### Human sCD146 ELISA Kit (orb406698)

**Description** Human soluble cluster of differentiation 146,sCD146 ELISA Kit

**Reactivity** Human

**Range** 1.25 ng/mL-80 ng/mL

**Target** sCD146

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [P43121](#)

**Sensitivity** 0.31 ng/mL

**Expiration Date** Please enquire.