

Product Datasheet

Human Galectin 3 ELISA Kit (orb406712)

Description	human lectin-galactose binding-soluble 3,Lgals3 ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	LGALS3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, saliva
Uniprot ID	P17931
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.