

Product Datasheet

Human ABCC5 ELISA Kit (orb406719)

Description	Human Soluble Mesothelin Related Peptide(SMRP)ELISA Kit
Reactivity	Human
Range	0.312 pmol/ml-20 pmol/ml
Target	ABCC5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O15440
Sensitivity	0.078 pmol/ml
Expiration Date	Please enquire.