

Product Datasheet

Human TNFRSF1B ELISA Kit (orb406725)

Description Human soluble tumor necrosis factor receptor 2,sTNF-R2 ELISA Kit

Reactivity Human

Range 31.25 pg/mL-2000 pg/mL

Target TNFRSF1B

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID **P20333**

Sensitivity 7.8 pg/mL

Expiration Date Please enquire.