

## Product Datasheet

### Human NOTCH3 ELISA Kit (orb406734)

<b>Description</b>	Human Neurogenic locus notch homolog protein 3(NOTCH3) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	125 pg/mL-8000 pg/mL
<b>Target</b>	NOTCH3
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>Q9UM47</b>
<b>Sensitivity</b>	31.25 pg/mL
<b>Expiration Date</b>	Please enquire.