

Product Datasheet

Human SERPINB3 ELISA Kit (orb406785)

Description	Human squamous cell carcinoma antigen 1 (SCCAg1) ELISA kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Target	SERPINB3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum,plasma,tissue homogenates.
Uniprot ID	P29508
Sensitivity	0.156 ng/ml
Expiration Date	Please enquire.