

## Product Datasheet

### Human Leptin Receptor ELISA Kit (orb406833)

<b>Description</b>	Human Leptin receptor,LR/Ob-R ELISA KIT
<b>Reactivity</b>	Human
<b>Range</b>	3.12 pg/mL-200 pg/mL
<b>Target</b>	LEPR
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P48357</a>
<b>Sensitivity</b>	0.78 pg/mL
<b>Expiration Date</b>	Please enquire.