

## Product Datasheet

### Human Annexin A3 ELISA Kit (orb406839)

<b>Description</b>	Human Annexin III(ANX-III) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.78 ng/mL-50 ng/mL
<b>Target</b>	ANXA3
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<a href="#">P12429</a>
<b>Sensitivity</b>	0.195 ng/mL
<b>Expiration Date</b>	Please enquire.