

## Product Datasheet

### Human OSCAR ELISA Kit (orb406904)

<b>Description</b>	Human Osteoclast-associated immunoglobulin-like receptor(OSCAR) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	6.25 ng/mL-400 ng/mL
<b>Target</b>	OSCAR
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q8IYS5</a>
<b>Sensitivity</b>	1.56 ng/mL
<b>Expiration Date</b>	Please enquire.