

## Product Datasheet

### Human PGRN ELISA Kit (orb406919)

<b>Description</b>	Human progranulin(PGRN) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	4.7 ng/mL-300 ng/mL
<b>Target</b>	PGRN
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Sensitivity</b>	1.17 ng/mL
<b>Expiration Date</b>	Please enquire.