

Product Datasheet

Human LTBP2 ELISA Kit (orb406933)

Description	Human latent transforming growth factor beta binding protein 2,LTBP2 ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	LTBP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q14767
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.