

## Product Datasheet

### Human CCL26 ELISA Kit (orb406950)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human C-C motif chemokine 26(CCL26) ELISA kit |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 39 pg/ml - 2500 pg/ml                         |
| <b>Target</b>            | CCL26   |
| <b>Note</b>              | For research use only                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates             |
| <b>Uniprot ID</b>        | <a href="#">Q9Y258</a>                        |
| <b>Sensitivity</b>       | 9.75 pg/ml                                    |
| <b>Expiration Date</b>   | Please enquire.                               |