

Product Datasheet

Human SIRT6 ELISA Kit (orb406961)

Description	Human NAD-dependent deacetylase sirtuin-6 (SIRT6/SIR2L6) ELISA kit
Reactivity	Human
Range	3.12 pg/mL-200 pg/mL
Target	SIRT6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q8N6T7
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.